




**PMH** PRODUTOS  
MÉDICOS  
HOSPITALARES



Standard

## **BLOOD ADMINISTRATION SYSTEM WITH 200M 25CM<sup>2</sup> FILTER, 160CM OF TUBE, INJECTION POINT AND FIXED LUER LOCK**

201.0250

### Technical description

Blood delivery system with 200? M 25cm<sup>2</sup> filter,  
160cm tube, injection point and fixed luer lock

REFERENCE	201.0250
CLASS	Class II a
STERILITY	ESTÉRIL
EAN	5607716002509

**Care with excellence since 1990**

*research and development medical products tailored  
to customer needs hospital*

## TECHNICAL FEATURES

Gravity blood delivery system with:

Cover for punch protection

Single-channel punch for easy punching in bags, placed in the chamber

Phthalate-free PVC chamber, transparent and flexible allowing continuous drip observation

200cm filter of 25cm<sup>2</sup>

PVC tube free of phthalates and DEHP, transparent, 3.00x4.10mm, 160cm

High-precision flow control ram, a design that allows you to maintain the permeable path, greater adaptability, greater safety, with a safe system for disposal.

Injection point for latex-free rubber needle

Luer lock male fixed

Protective cover luer lock

## ADDITIONAL INFORMATION

### COMPOSITION

ABS, LDPE, PVC, PA, POLISOPRENE  
FREE FROM DEHP AND PHALATES  
LATEX FREE

### BOX CHARACTERISTICS

### CONSERVATION

TEMPERATURE BETWEEN 0 ° C AND + 40 ° C

### CHARACTERISTICS OF THE PACKAGING

UNIT PACKAGING IN PEEL PACK, STERILE-ETO, EASY TO OPEN, WITH INDICATION OF THE PRODUCT REFERENCE, MANUFACTURING DATE, BATCH NUMBER, EXPIRY DATE, MARKING CE0086, INSTRUCTIONS FOR USE AND APPLICABLE SYMBOLOGY

### WARNINGS AND PRECAUTIONS


STERILE SINCE NOT OPENED OR DAMAGED. APIROGENIC.  
USE ONLY IF THE PROTECTIVE COVERS ARE PLACED

### DIRECTIONS FOR USE

ADMINISTERING BLOOD TO A PATIENT FROM A BLOOD BAG

### EXPIRATION DATE

5 Years



[WWW.PMH.PT](http://WWW.PMH.PT) | TL. +351 263 650 680

PMH, SA - PRODUTOS MÉDICO HOSPITALARES  
ZONA INDUSTRIAL DA MURTEIRA, RUA GUINÉ BISSAU, LOTE 9  
2135-327 SAMORA CORREIA